

Materials and Equipment

Materials

- Cardiac Progenitor Cells ([AC015-F](#))
- Cardiac Cellutions Medium ([C-MGRO](#))
- Cardiomyocyte Cellutions Differentiation Medium ([C-MDIFF](#))
- FGF-2 (20ng/ml, concentration in C-MGRO, to be added – not supplied)
- 5-azacytidine (5-AZA – not supplied; to add when specified: 10 µM in C-MDIFF)
- TGF-β-1 (not supplied; to add only when specified: 1 ng/ml in C-MDIFF)
- Collagen Coating Solution (150 µg/ml)
- 1× PBS (Ca²⁺/Mg²⁺-free)
- Fetal Bovine Serum (FBS)
- TrypLE Express or equivalent enzyme solution
- 15.0 ml and 50.0 ml conical centrifuge tubes
- Cell culture vessels
- Dimethyl Sulfoxide (DMSO)

Equipment

- Biosafety Cabinet
- Water Bath, 37°C
- Centrifuge
- Incubator, 37°C, 5% CO₂
- Pipet-Aid
- Pipettes

Procedure

1.	<p>Immediately upon delivery</p> <ul style="list-style-type: none"> • Remove vial from shipping container and check that it is still frozen—if thawed, contact DV Biologics. • Transfer frozen vial immediately to liquid nitrogen until you are ready to thaw and begin cell culture. • Prior to thawing cells, make sure to have all required reagents and materials ready to use. • Coat cell culture vessels with 150 µg/ml collagen coating solution (0.2 ml/cm²). <ul style="list-style-type: none"> ○ Incubate at room temperature for 2 hrs. ○ Aspirate collagen coating solution. ○ Wash the culture vessel with PBS.
2.	<p>Thawing Cells</p> <ul style="list-style-type: none"> • Perform and maintain all cell culture using aseptic techniques. • Aliquot and warm only media required. To the C-MGRO for use that day, add fresh FGF-2 (20ng/ml). • Thaw cells rapidly and with agitation (See How to Thaw Cells). <ul style="list-style-type: none"> ○ Slowly dilute thawed cells with 3- to 10-fold excess cold FBS (or complete C-MGRO medium). ○ Centrifuge at 400 g for 5 minutes. • Remove supernatant. Disperse cell pellet in 1-5 ml fresh 37°C warm C-MGRO medium with FGF-2. • Count Cells (See How to Count Cells).
3.	<p>Plating and Maintaining Cells</p> <ul style="list-style-type: none"> • Carefully add cell suspension to collagen-coated culture vessel. Seed at density of 1.0x10⁴ cells/cm². • Leave cells undisturbed for 2 days. Feed cells on the third day by completely replacing the medium. NOTE: There may be many unattached cells. These cells could be progenitors or stem cells that have not yet attached. You may collect these, centrifuge and reseed. • Feed cells every 2-3 days with fresh C-MGRO medium (with FGF-2 added) by either: <ul style="list-style-type: none"> ○ Performing a complete media replacement and discarding the spent media, or ○ Removing the spent media but not discarding until after collecting any unattached cells by centrifugation as noted above. After centrifuging, aspirate supernatant and disperse pellet in C-MGRO with FGF-2 to reseed. • Once cells become 60-70% confluent, feed more regularly (every other day). • When cells are nearing confluency (approximately 80%), prepare to passage cells either for propagation/expansion of the progenitor cells or for their terminal differentiation. NOTE: Do not let cells become overly confluent.

4.	<p>Passaging of Cells (for Details see How to Passage Cells)</p> <ul style="list-style-type: none"> Aspirate and discard medium; wash cells with PBS without Ca²⁺ and Mg²⁺, and aspirate. Add TrypLE Express or equivalent enzyme solution at a volume that just covers the cell monolayer. Incubate at 37°C; microscopically observe detachment of cells every 5 minutes; avoid toxic effects of enzymatic over-digestion. Quench the enzyme digestion by adding 2- to 5-fold volume of C-MGRO medium. Transfer the entire resulting suspension of detached cells from the culture vessel to a centrifuge tube. Centrifuge the cell suspension at 400 g at 4°C for 5 min. Aspirate supernatant and disperse cell pellet in warm C-MGRO medium for counting. Count cells (see How to Count Fresh Cells). To propagate cardiac progenitor cells, seed cells at 1.0 x 10⁴ cells/cm² following the differential adhesion/pre-plate method described in the next section. If more cells are obtained than needed for re-plating, additional cells may be frozen. (See How to Freeze Cells).
5.	<p>Differential Adhesion/Pre-Plate Method</p> <ul style="list-style-type: none"> Coat flask with collagen and set aside until time for second seeding of cells. For the first seeding of cells, plate cells on NON-coated flasks at seeding density of 1.0 x 10⁴ cells/cm² in C-MGRO medium. Incubate for 30 minutes in 37°, 5% CO² incubator. Transfer medium and non-adhered cells to a tube. Rinse flask with fresh C-MGRO medium and pool rinse with non-adhered cells. Centrifuge at 400 g for 5 minutes at 4°C. Aspirate supernatant and disperse cell pellet in C-MGRO medium for counting. Count cells (see How to Count Fresh Cells). Seed cells into collagen-coated flask at a density of 1x10³ to 5x10³ cells/cm². The non-coated flask, which has a mix of cells including some fibroblasts and cardiac muscle cells, can be discarded. Repeat these steps for each passage until ready to differentiate. For terminal differentiation use p3 to p5 cells.
6.	<p>Terminal Differentiation</p> <ul style="list-style-type: none"> Seed cells in flask/wells using complete growth medium: C-MGRO (with FGF-2 added), prepared as described previously. Once confluent, feed cells with complete differentiation medium: C-MDIFF, consisting of basal combined with supplement – without any other added component at this step. NOTE: It is important that cells be precisely 100% confluent. Exceeding or not fully reaching 100% confluency could affect success of differentiation. Let cells acclimate to new C-MDIFF medium for 1 day. (Do not use FGF-2 or TGF-β-1 for this step). This will be Day 0. On Day 1, feed cells with a complete medium replacement using C-MDIFF with 5-AZA. Prepare this by warming fresh aliquot of C-MDIFF and adding 5-AZA (10 μM, final concentration). Repeat the above step for Day 2 and Day 3, feeding cells by replacing with fresh C-MDIFF with 5-AZA (10 μM) added fresh each time. On the fourth day, prepare fresh C-MDIFF with TGF-β-1 WITHOUT 5-AZA. Prepare this by warming fresh aliquot of C-MDIFF and adding TGF-β-1 (1ng/ml, final concentration). Maintain cells in C-MDIFF with TGF-β-1 from Day 4 onward for 2-3 weeks, feeding cells every 2-3 days with fresh C-MDIFF with TGF-β-1. NOTE: Observe cells regularly for lifting and toxicity.