

Materials and Equipment

Materials

- Cardiomyocytes ([AC008-F](#))
- Cardiomyocyte Cellutions Maintenance Medium ([C-MAIN](#))
- TGF- β -1 (1 ng/ml, concentration in C-MAIN, to be added – not supplied)
- Collagen coating solution (150 μ g/ml)
- 1 \times PBS (Ca^{2+} / Mg^{2+} free)
- Fetal Bovine Serum (FBS)
- 15 ml, 50 ml conical centrifuge tubes
- Cell culture vessels

Equipment

- Biosafety Cabinet
- Water Bath 37°C
- Centrifuge
- Incubator, 37°C, 5% CO₂
- Pipet-Aid
- Pipettes

Procedure

1.	<p>Immediately Upon Delivery</p> <ul style="list-style-type: none"> • Remove vial from shipping container and check that it is still frozen—if thawed, contact DV Biologics. • Transfer frozen vial immediately to liquid nitrogen until you are ready to thaw and begin cell culture. • Prior to thawing cells, make sure to have all required reagents and materials ready to use. • Coat cell culture vessels with 150 μg/ml collagen coating solution (0.2 ml/cm²). <ul style="list-style-type: none"> ○ Incubate at room temperature for 2hrs. ○ Aspirate collagen coating solution. ○ Wash the culture vessel with PBS.
2.	<p>Thawing Cells</p> <ul style="list-style-type: none"> • Perform and maintain all cell culture using aseptic techniques. • Prepare C-MAIN medium. After mixing, store complete medium at 4°C for a maximum of 30 days. • Aliquot and warm only media required for use that day at 37°C. • To the C-MAIN medium for use that day, add fresh TGF-β-1 (1 ng/mL). • Thaw cells rapidly and with agitation (See How to Thaw Cells). <ul style="list-style-type: none"> ○ Slowly dilute thawed cells with 3- to 10-fold excess cold FBS (or complete C-MAIN medium). ○ Centrifuge at 400 g for 5 minutes. • Remove supernatant. Disperse cell pellet in 1-5 ml fresh 37°C warm C-MAIN medium with TGF-β-1. • Count Cells (See How to Count Cells).
3.	<p>Plating and Maintaining Cells</p> <ul style="list-style-type: none"> • Carefully add cell suspension to collagen-coated culture vessel. <ul style="list-style-type: none"> ○ Plate cells at a minimum cell seeding density of 5.0 x 10⁴ cells/cm². ○ Leave cells undisturbed for a minimum of 2 days. • Feed cells on the third day by completely replacing the medium with fresh C-MAIN with TGF-β-1. NOTE: There may be many unattached cells. You may collect these, centrifuge and reseed. • Maintain cells in C-MAIN with TGF-β-1, feeding every 2-3 days, and observe cells daily. <ul style="list-style-type: none"> ○ Cells should be fed and maintained in culture for a minimum of 10 days prior to use. ○ Cells should be kept confluent for optimal results.